



RADIALPET

Radial Wave Therapy (RWT) for Pet Orthopedic Disorders

Radialpet, the ultimate solution for superficial chronic orthopedic indications, is a Radial Wave Therapy (RWT) System. Radialpet was developed after extensive theoretical knowledge and vast practical experience in the field of Extracorporeal Shock Wave Therapy (ESWT) and Radial Wave Therapy (RWT), and is considered one of the leading devices in the Orthopedic field.

RWT has demonstrated to be a valuable treatment for a variety of musculoskeletal injuries and manage to treat superficial pathologies where conventional methods have failed, with an extremely high success rate as ESWT.

When other therapies, such as corticosteroids and conventional physiotherapy, fail, Radialpet delivers highly efficient results. It is a safe, non-invasive and rapid non-surgical method for treating resistant or non-responsive superficial orthopedic pathologies, improving pain and function level dramatically.

RWT - up to 85% success rate on chronic indications!

- Extremely low maintenance costs
- No anesthesia required
- Comfortable treatment
- Short treatment time
- Flexibility and ease of use
- Compact and mobile
- High localization precision with no imaging required



GROUNDBREAKING TECHNOLOGY Within Your Reach



How Radialpet works?

The radial waves are applied directly onto the pain source. Radial waves act physically on the tissue and induce multiple favorable biological effects:

- Attenuate inflammatory response or attenuates inflammation
- Vasodilation, increase in blood flow and oxygen supply to the damaged tissue
- Release of ANGIOGENIC factors leading to angiogenesis and neovascularization
- Tissue healing
- Fragmentation of calcium deposits into “crumbs” subsequently absorbed into the tissue
- Differentiation of mesenchymal stem cells into bone forming cells (Osteoblasts) which leads to bone restoration and fracture healing
- Analgesic effect

The treatment protocol varies between four to six sessions in average, depending on the animal’s condition, the severity of the pathology and the individual healing rate.

An outpatient procedure lasts up to 5-10 minutes.

Radialpet unit includes three different treatment tips to optimize pain relief and ensure a precise treatment area.

Multiple applications

- Joint disease, including:
 - Hip arthritis (coxarthrosis)
 - Knee arthritis (gonarthrosis)
 - Elbow arthritis (cubarthrosis)
- Osteoarthritis, including:
 - Hip dysplasia, hip osteoarthritis
 - Elbow osteoarthritis
 - Knee osteoarthritis
 - Shoulder osteoarthritis
- Tendinopathies, including Tendinopathy of biceps tendon



203 Perry Parkway, Suite #6, Gaithersburg, MD 20877, USA
Tel: +1 (301) 944-1575
marketing@medispec.com
www.medispec.com

