



4-Axis Urology Treatment Table

| Feature | 4-Axis Treatment Table | |
|---------------------|---|---------------|
| Movements | Height: | 7.8" |
| | X Axis (head-to-toe): | 3.9" |
| | Y Axis (side-to-side): | 3.9" |
| | Trendelenburg | 7° |
| | Reverse Trendelenburg: | 7° |
| | Lateral Roll: | None |
| Dimensions | Height: | 34.6" – 42.5" |
| | Width (tabletop): | 20.5" |
| | Length (fully extended): | 87.8" |
| Max Load | 400 pounds distributed load | |
| Construction | <ul style="list-style-type: none"> • Carbon fiber table top with dual cut-outs • 2" mattress • 4 segments of accessory rails • Detachable leg and head extensions • Front Wheels: 3" Swivel casters • Rear Wheels: 4" Swivel casters with locks | |
| Controls | <ul style="list-style-type: none"> • Hand Control – all motions | |
| Power | <ul style="list-style-type: none"> • 120 V / 60 Hz single phase | |

*Specifications are subject to change without notice.



20410 Observation Drive, Suite # 102, Germantown, MD 20876-6419
P: (301) 944-1575 . Toll Free: (888) 663-3477 F: (301) 972-6098
e-mail: info@medispec.com www.medispec.com