



5-Axis Endourology Treatment Table

Feature	5-Axis Treatment Table	
Movements	Height:	10"
	X Axis (head-to-toe):	4"
	Y Axis (side-to-side):	10"
	Trendelenburg	0° - 20°
	Reverse Trendelenburg:	20°
	Lateral Roll:	20°
Dimensions	Height:	33.5" – 43.5"
	Width (tabletop):	24"
	Length (fully extended):	93"
Max Load	450 pounds distributed load	
Construction	<ul style="list-style-type: none"> • Carbon fiber table top with dual cut-outs • 2" mattress • 4 segments of accessory rails • Detachable leg extension • Front Wheels: 3" swivel casters with locking system • Rear Wheels: 5" swivel casters with central locking 	
Controls	<ul style="list-style-type: none"> • Hand Control – all motions • Foot Control – Trendelenburg & lateral 	
Power	<ul style="list-style-type: none"> • 120 V / 60 Hz single phase • Optional battery back-up 	

*Specifications are subject to change without notice.