





ENDO COMFORT™ TREATMENT TABLE #1 IN VERSATILITY

Suitable for a
Range of Medical
Procedures



Mobile



Stationary

STANDALONE UNIT

DESIGNED FOR FLEXIBILITY AND COMFORT

One Versatile Table for Different Functions
Compatible with most portable lithotripters, x-rays, and C-arms

Endo Comfort™ provides a flexible, comfortable treatment table, available as a mobile or stationary unit. It accommodates a variety of procedures, including endo-urology, lithotripsy, radiology, and a range of other treatments.




For a Variety of Endo-urology & Lithotripsy Procedures
Customizations to meet patient and doctor needs


Endo Comfort™ is an ideal table for most any clinic, hospital, or medical offices, Endo Comfort™'s advantages include its small footprint, convenient mobility, and easy access to the patient on both sides of the treatment table.



ZERO CALIBRATION



ZERO SETUP TIME



ZERO FOOTPRINT

EASY PATIENT POSITIONING

MULTI-DIRECTIONAL MOVEMENTS FOR EASY WORKFLOW

- Using the C-arm, high resolution images can be taken from any angle
- All movements are controlled via a remote control unit
- Voltage (AC) Single Phase 100-240 volts



5 DIRECTIONS

3-AXIS, 4-AXIS, AND 5-AXIS OPTIONS




- Longitudinal
- Lateral
- Vertical
- Endo Tilt (Trendelenburg)
- Endo Rotate

ENABLE CLEAR DIAGNOSIS

CARBON FIBER FOR RADIOTRANSPARENCY

- The transparent carbon fiber material allows seamless workflow
- X-ray images can be seen and captured in real time, without obstructions or limitations
- 4-part PVC customizable, flat mattress: head, main, legs, notch cover
- The main part contains left and right notches to enable treatment access from either side



ENDO COMFORT™



**Detachable Legs/Head
Extensions**



Mount Rail



Drain Bag

COMPREHENSIVE ACCESSORIES

ADJUST IN REAL-TIME

- Flat mattress with four parts and a head pillow
- Lateral mount rails for stirrups, IV poles, patient restraints, and other equipment
- Supporting plates for notch covers
- Detachable leg & head extensions
- Hand rail
- Arm rests (optional)
- Infusion fluid holder (optional)

Since its founding 30 years ago, Medispec has been a leader in the advancement of wave therapy technologies, initiating new treatments and applications across medical disciplines, providing health and economic benefits to both providers and their patients.



LIVE BETTER. LOVE BETTER.



203 Perry Parkway, Suite 6,
Gaithersburg, MD 20877, USA
Tel: +1 (301) 944-1575
Fax: +1 (301) 972-6098
E-mail: marketing@medispec.com
www.medispec.com