

# ENDO-UROLOGY TABLE™

## 5-axis Multipurpose Motorized Treatment Table

Medispec's Endo-Urology table is your best choice of treatment Table designated for ESWL as well as other Endo-Urology treatments.

Medispec's 5-axis table is designed to optimize a variety of endourological procedures as well as other specialties by allowing complex procedures to be performed.

The Treatment Table has cutout on both sides to accommodate the Lithotripter shockwave head.

The tabletop is made of carbon fiber for obstruction-free imaging.

Medispec's mobile treatment table has a weight capacity of 400 lbs (181Kg).




- Compatible for endourologic suites
- Comfortable and mobile
- Effortless positioning
- Very low maintenance
- Transparent for X-ray



## The 5 Axes Motion

The medispec's treatment table has a motor-driven movement for Trendelenburg/rotation tilt and horizontal/vertical motion.

Positioning the patient is effortless with ergonomic hand held remote control.



## Treatment Table Accessories

A selection of accessories maximizes the treatment table usability for both ESWL and various endourologic procedures:

- 4-parts flat mattress and a head pillow
- Supporting plate for the notch cover
- Detachable legs/head extensions
- Hand rail
- Arm rests (optional)
- Detachable legs stirrups (optional)
- Infusion holder (optional)

